Thera **Lite**

*Aura*

# Mood and Energy Enhancing Light User Guide

---

**⚠ WARNING:** This light should not be used on patients with macular degeneration. **ALWAYS** consult a physician before using this product.

---

**COMPASSHEALTH**

Manufactured for: Compass Health Brands Corp.  
Middleburg Heights, OH 44130  
[www.compasshealthbrands.com](http://www.compasshealthbrands.com) • Phone 800.376.7263  
Made in China • 42-FGP83000 0000 02 © 2021

Thank you for purchasing the TheraLite Aura. Please read this manual in its entirety. For more information on light therapy, please visit [carex.com](http://carex.com).

TheraLite Aura has been designed with your needs in mind. Projecting light from above, adjustable angle and multiple light settings make the TheraLite Aura the perfect companion at home or at the office.

UV-free LED technology and a comfortable color temperature ensure you experience a safe, comfortable therapeutic light. The adjustable lamphead lets you switch from therapy lamp to desk light with the turn of the knob.

### **Bright Light Benefits**

Daily exposure to bright light helps to maintain your circadian rhythm, boost energy and manage sleeping patterns. Light therapy is also the first line of recommended treatment for the Winter Blues.

For many, morning sun exposure is limited by work, school or other commitments, even in the summer. The TheraLite is a bright light therapy lamp designed to fill this gap.

### **Included Items**

TheraLite Aura lamp, Philips screwdriver and User Guide.

### **Safety Precautions**

1. If you have eye problems or have had eye surgery, you should consult your doctor or ophthalmologist before you use any form of light therapy. Certain eye and eye-related conditions may respond unfavorably to light therapy.
2. As with any electrical device, **DO NOT** use the TheraLite in or near water.
3. **DO NOT** overuse the TheraLite. Overuse may cause irritability, excessive energy, or difficulty falling asleep at bedtime.
4. It is **NOT** necessary to stare into the light. Read, work or talk on the phone during your light therapy session.
5. During the first few uses, your TheraLite may give off an odor; this is harmless, will not affect the use and will quickly disappear.
6. **ALWAYS** use in a well-ventilated area away from other heat sources and ensure that air can flow freely around the lamphead when the product is in use.

### **Disclaimer**

TheraLite is an innovative light supply system and is not a listed medical device in the USA.

### **Warranty**

Your lamp is maintenance-free and should provide years of service. The TheraLite has a one-year limited warranty from the date of purchase by the original purchaser. Please retain proof of purchase to ensure eligibility for warranty coverage.

Compass Health Brands, may choose to repair or replace a defective unit. Shipping costs are the responsibility of the product owner.

Misuse, abuse, accidental or incorrect operation, contrary to the instructions, will void a warranty claim.


## Specifications

Power:	120V, 60Hz, 35W
Size:	22" x 11.25" x 11.5" / 56cm x 28cm x 29cm
Lens:	High-impact polycarbonate
Light Intensity:	10,000 LUX at 12"/30cm
Light Source:	168 LED Bulbs, 3000K, 100% UV Free
Angles:	Adjustable

## Light Therapy Basics


LUX or the light intensity at a specific distance from the light is the common term used in bright light therapy. The chart below details the average treatment times for the Winter Blues based on your eyes proximity to the lamp.

DISTANCE	SESSION TIME
12"/30cm	30 minutes
18"/45cm	60 minutes
24"/60cm	120 minutes

## Assembly Instructions

1. On a level, solid surface lay out the lamphead/arm, the lamp base and the 2 Philips screws for assembly.
2. Position the lamphead/arm on its side.
3. Carefully and firmly insert the bottom of the lamphead/arm into recess on the lamp base in the orientation shown.
4. Insert and by hand seat the Philips screws through the holes on the bottom of the lamp base.
5. Using the Philips screwdriver, hand-tighten the screws until snug – be careful not to overtighten.
6. Stand the lamp upright.
7. While holding the lamphead with one hand, with the other, turn the angle adjustment knob to loosen the lamphead, position to the desired angle and re-tighten the adjustment knob.
8. Plug the lamp into a standard 120V electrical outlet
9. Use the on/off touch button to turn the lamp on. Press "+" to increase the brightness and press "-" to lower the brightness.

### Using Your TheraLite for Therapy


For those with mild to moderate Winter Blues, sleep or circadian rhythm issues, regular daily use of the TheraLite can help provide relief. For most conditions, use shortly after waking or in the morning.

Place the TheraLite Aura above your line of sight, where possible, to maximize the light reaching your eyes. You should be facing the lamp with eyes open during sessions, but you do not need to look directly at the light.

Over the first few days of use, assess the effect and, if necessary, lengthen or shorten the session time based on your response.

### Symptoms of Overuse

Negative effects are rare with light therapy. However, if you experience increased irritability, excessive energy, and/or any consistent discomfort during or after use, decrease session time or move the light further away. Using the lamp late in the afternoon or in the evening is not advised, as this may disrupt your sleep. It is recommended that you consult your doctor before beginning any bright light regimen.

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Fabriqué pour : Compass Health Brands Corp.