

# Invacare® Top End® Force CC™ Handcycle

## Price List and Order Form

Effective April 15, 2012

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.topendwheelchair.com



Provider Name \_\_\_\_\_ Date \_\_\_\_\_ Top End Job Number \_\_\_\_\_  
 Account # \_\_\_\_\_ P.O. # \_\_\_\_\_ Phone # \_\_\_\_\_ Fax # \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_  
 Customer Name \_\_\_\_\_ Phone No \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_  
**Ship order to:**  Provider  Customer **Order will ship to Provider unless otherwise noted.**

### TOP END FORCE CC HANDCYCLE:

FRCCC Top End Force CC Off Road 27 Speed Handcycle \$4900.00

#### Seat Width: (Select One)

- 15" Stealth Black Frame and Fork  
 17" Stealth Black Frame and Fork

#### Braking System: Front Disk Brake and Parking Brake

XLT88 Dual Rear Disk Brake upgrade .....500.00

#### WHEELS: 559-26" wheels w/threaded Axles (Select Tire)


- Knobby mountain bike tires.....Standard  
 Cruiser clincher tires..... No charge

### SAFETY & COMFORT OPTIONS

NOTE: All non-standard options will be shipped separately

- |                                  |  |          |
|----------------------------------|--|----------|
|                                  | Seat and Back Cushions.....  | Standard |
|                                  | Chain Guard.....   | Standard |
|                                  | Safety Flag.....   | Standard |
| <input type="checkbox"/> XLT9    | Leg Guard Attachment.....  | 75.00    |
| <input type="checkbox"/> XLT84   | Crank width adaptor set (adds 1.5" per set) .....  | 20.00    |
| <input type="checkbox"/> XLT71   | Tool and Tire Repair Kit.....  | 90.00    |
| <input type="checkbox"/> XLT12   | Mirror, Leg Rest Mount.....  | 30.00    |
| <input type="checkbox"/> XLT13   | Water Bottle & Cage.....   | 35.00    |
| <input type="checkbox"/> XLT28   | Backpack Hydration System.....   | 100.00   |
| <input type="checkbox"/> XLT141  | Heart Rate Monitor.....  | 100.00   |
| <input type="checkbox"/> XLT14   | Cordless computer.....   | 125.00   |
| <input type="checkbox"/> XLT15   | Helmet.....  | 75.00    |
| <input type="checkbox"/>         | Small (6 5/8" - 7") <input type="checkbox"/> Medium (6 7/8" - 7 1/2") <input type="checkbox"/> Large (7 1/4" - 7 5/8") |          |
| <input type="checkbox"/> XLT16   | Extra Safety Flag.....   | 30.00    |
| <input type="checkbox"/> XLT6    | Front Safety Light.....  | 25.00    |
| <input type="checkbox"/> XLT7    | Rear Flashing Light.....   | 25.00    |
| <input type="checkbox"/> XLT 33  | Alignment Gauge.....   | 75.00    |
| <input type="checkbox"/> XLT20   | Bike Rack for Trailer Hitch.....   | 450.00   |
| <input type="checkbox"/>         | <input type="checkbox"/> 1 1/4" (Std.) <input type="checkbox"/> 2"   |          |
| <input type="checkbox"/> XLT40   | Cable Kit.....   | 35.00    |
| <input type="checkbox"/> XLT74   | Extra Steering Dampener.....   | 10.00    |
| <input type="checkbox"/> XLT86   | LeMond® Revolution™ Trainer w/ power meter.....  | 1250.00* |
| <input type="checkbox"/> XLT87   | LeMond Revolution Trainer w/o power meter.....   | 850.00*  |
| <input type="checkbox"/> XLT42   | Indoor Training Roller.....  | 250.00   |
| <input type="checkbox"/> TER126  | Click Strap Large (pad length~12").....  | 86.00    |
| <input type="checkbox"/> TER125  | Click Strap Med (pad length~10").....  | 80.00    |
| <input type="checkbox"/> TER124  | Click Strap Small (pad length~8").....   | 75.00    |
| <input type="checkbox"/> RACE22* | Top End T-Shirt <input type="checkbox"/> L <input type="checkbox"/> XL.....  | 10.00    |
| <input type="checkbox"/> RACE4*  | Wheel Bag.....   | 180.00   |
| <input type="checkbox"/> RACE20* | Handcycle Travel Bag.....  | 300.00   |
| <input type="checkbox"/> RACE13* | Top End Bike Jersey <input type="checkbox"/> L <input type="checkbox"/> XL.....  | 70.00    |

\* Standard Dealer discount does not apply, additional shipping charges apply.



### FORCE CC HANDCYCLE SPECIFICATIONS

**Frame** – Ovalized 7000 series aluminum, Stealth Black frame color

**Wheel base:** 60" Standard, Overall Length 86"

**Removable Fork** – 7000 series aluminum with adjustable crank housing.

**Components:** Shimano® Deore Rapid-fire® Brake/Shifter

**Mounting Position: Shifter / Brake Lever:** Right Handle

**Front Wheel Cassette:** 11-34

**Chain ring:** 44-34-22

**Handpedal Size:** Medium (1 3/4")

**Crank arms:** 170mm length x 14" wide

**Diagonal and Vertical Crank Height Adjustment**

**Camber: 0 degrees**

**Turning Radius** – 18 feet for 360° circle, 12 feet w/one backing maneuver.

**Static Tip Angle** - 39° for 250 lb subject.

**Transfer Seat Height** – 11.5" with 2" cushion.

**Ground clearance-** 6.5"

**Backrest** – Adjustable fore and aft.

**Back Angle** – Adjustable angle back of 90° - 120°.

**Propulsion System** – Shimano components 27 speeds, 44-34-22 chain rings and 11/34 cassette. The 11/34 cassette is operated by the Hands-On Rapid Fire unit which is mounted on the right hand pedal. It operates both the 11/34 cassette and the hands-on front disk brake. The top chain rings utilizes a cable driven shifter system located on upper crank. A parking brake/second brake is mounted on the upper fork which also operates the optional rear disk brakes if ordered.

All prices subject to change w/o notice. All prices in U.S. Dollars. All orders F.O.B. Pinellas Park, Florida, USA

For the most current pricing information go to www.topendwheelchair.com.

# Invacare® Top End® Force CC™ Handcycle

**Crank set** – Top End Crankset

**Steering** – Pivot steering and raked fork for controlled turning and stability.

**Crank** – Adjustable height both vertical and diagonal.

**Hand pedals** – Hand pedals are made of lightweight aluminum and are ovalized to precisely fit the hand. The bearings allow the pedal to rotate while pedaling are located in the crank set and not in the pedal, which makes this pedal very comfortable.

**Footrest** – Multi-adjustable fore and aft w/straps to securely cradle calves and feet. Straps are 4" wide and feature safety straps for ankles.

**Performance wheels** –Silver spoked wheels with threaded axles and alloy rims. 14 gauge stainless steel spokes. 559 (26)" x 1.75" Knobby Clincher Tire/Tube 40-60psi. or optional Cruiser Clincher Tire.

**Upholstery** – Black nylon upholstery with 2" seat and back cushion standard. Seat restraint included.

**Warranty** –3 year on frame and fork to original owner/limited parts warranty.

**Weight Limit:** 250 pounds

**Shipping weight and freight charges:** Shipping weight is 62 pounds. Box size is 77" x 27" x 23". Freight charge to ship within the continental USA is \$100 to a business address and \$155 to a residential address. Extra freight charges for large options such as the wheel bag, handcycle travel bag, training roller or bike rack may apply.



**Yes, you can.®**

**Invacare Corporation** [www.invacare.com](http://www.invacare.com)

**TOP END**  
Sports and Recreation  
Products  
4501 63<sup>rd</sup> Circle North  
Pinellas Park, FL USA  
33781-5914

---

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are registered trademarks of Shimano Inc. LeMond and Revolution are a registered trademarks of Lemond Fitness, Inc.  
©2012 Invacare Corporation

Form No. 12-291 Rev 04.12  
P2012-291