

BAJA® WRANGLER® 2




Shown in Black Camo

PROUDLY MADE IN
AMERICA
OF FOREIGN AND DOMESTIC PARTS

FDA Class II Medical Device**

THE SPIRIT OF ADVENTURE

This rugged outdoor explorer provides you with state-of-the-art suspension and drivability to make you feel like you're floating on air, industry leading dual-drive off-road capabilities with independent traction and torque the competition will wish they had. The Baja® Wrangler® 2 is packed with the go-power of up to 11.4 mph and built to get you to the places you want to go and the places you haven't thought of yet.



BAJA® MOBILITY

BAJA® WRANGLER® 2

FEATURES

K0899 POWER MOBILITY DEVICE:

Not coded by DME PDAC or does not meet criteria



FDA Class II Medical Device**

Shown with optional front basket

- CTS Suspension includes industry-leading front and rear suspension for a smoother ride
- Steering dampening system absorbs uncontrolled movement over rugged terrain, limits fatigue and reduces energy needed to control handles
- Choice of 16"x16-18" or 18"x18-20" black vinyl, limited recline memory foam seat with manual adjustment, heavy-duty seat post, sliders and height, width, depth, and angle adjustable armrests (optional seat sizes of 20"x18-20" and 22"x20-22" available)
- User friendly LED console with variable speed control displays time, temperature, miles driven, and trips taken plus horn and light buttons
- Adjustable delta tiller with ergonomic, wraparound handles, XLR charging port and USB charger for charging portable devices
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Full LED lighting package includes headlights, hazard lights, and directional signals
- Rear backup sensor LED lights indicate when you are too close to an object when in reverse
- Hydraulic hand brake with standard regenerative and electro-mechanical braking system
- Choice of five color options including Black Camo and White Camo.



User-Friendly Console

NEW
Storage pod accessory attaches to rear of seat.



ACCESSORIES

- 3" single or double cane/crutch holder
- 50/60/70" lap belt
- Cup holder
- Cell phone holder
- Forearm crutch holder
- Front basket
- Oxygen holder
- Quad cane holder
- Rear basket
- Saddlebag
- Storage pod
- Walker holder
- Wishbone crutch holder

SPECIFICATIONS

MODEL NUMBER: BA401

WEIGHT CAPACITY: 450 lbs. TURNING RADIUS³: 94.75"

MAXIMUM SPEED¹: Up to 11 mph at 450 lbs.
Up to 11.4 mph at 200 lbs.

DIMENSIONS³: Width: 32.5" Length: 61.25"

BATTERY OPTIONS: 75 AH or 100 AH (optional)

RANGE PER CHARGE^{1,2}: 75 AH: Up to 13.2 miles at 450 lbs.
75 AH: Up to 20.2 miles at 200 lbs.
100 AH: Up to 17.4 miles at 450 lbs.
100 AH: Up to 26.6 miles at 200 lbs.

TOTAL WEIGHT⁴: 274 lbs. (with standard seat and batteries)

WARRANTY

FRAME: Lifetime limited

DRIVETRAIN: 2-year limited

ELECTRONICS: 2-year limited

BATTERY: 6-months

HCPCS Code: K0899

** Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. Battery weight may vary based on battery manufacturer.

AUTHORIZED PROVIDER



bajamobility.com

