



## **EW-46 Specifications**

Drive up to 13 mph for up to 35 miles on a single charge with EWheels highly anticipated EW-46. This fun-to-drive electric scooter mobility scooter features front and rear suspension for a smoother ride, while 16" pneumatic tires allow you to drive effortlessly over varied terrain.

The EW-46 fully adjustable large captain's seat is one of the most comfortable seats around and it holds up to 400lbs!

An easy to read LED console shows estimated speed, mileage and battery display for comfortable navigation.

A full LED lighting package delivers exceptional visibility with directional turn signals, brake lights and ultra bright headlight with both rear and side LED running lights.

Dual rear-view mirrors allow you to see what's behind you, while the Electric Brake Assist hand-braking system ensures safe braking.

Also the sturdy front steel bumper protects you from the unexpected.

Speed:  
13mph maximum speed

Range:  
Up to 35 miles per charge

Weight Capacity:  
400 lbs.

Brakes:  
Electric Brake Assist

Tubeless Tires:  
Yes

Executive Seat:  
Fully Adjustable Seat,  
Rear Tilt & Adjustable Tilt Up Armrest

Seat Dimensions:  
18"x18"x28" (LxWxH)

Electronic Package:  
Digital Dash, Front &  
Rear Lights, Brake Lights, Turn Signal's,  
LED Rear Light Strip & Side Reflector  
Running Lights

Suspension:  
Front and Rear Spring  
Suspension

Bumper:  
Heavy Duty Metal Bumper

Rear Basket:  
Basket Size 10"x20"x11"  
(LxWxH)

Steering Height:  
42" From Floor

Turning Radius:  
73"

Motor:  
500 Watt Brushless Motor

Batteries:  
48 Volt, 20 AH

Wheelbase:  
40"

Dimensions:  
63"x28"x42" (LxWxH)

Floor Clearance:  
4.5"

High/Low Speed Switch:  
3 Setting Adjustable Speed Control Switch

Tire Size:  
Front and Rear 16" Tires

Charger:  
48 Volt Smart Charger

Throttle Type:  
Twist Throttle

Key Start:  
Yes

Scooter Weight:  
260 lbs

Battery Indicator:  
Yes, Digital

Warranty:  
Three Year Limited