



## Model EW-72 Spec Sheet

Scooter Weight Capacity  
500Lbs.

Speed:

Up to 15mph Distance: 43 miles per charge (distance varies on terrain, riders weight, road surface, etc.)

Front Seat:

Type Adjustable Folding Captains chair with backrest

Front seat:

Size 18 ½ "x 18 ½ "x 30" (LxWxH)

Front Seat:

26 ½ "from Ground

Carton Size:

Shipped Fully assembled on Pallet – 63"x 36"x 54" (LxWxH)

Basket Size:

10"x20"x12" (LxWxH)

Steering Height

46" from Ground

Turning Radius  
92"

Power:  
Electric

Watts:  
700 Watt

Motor Type:  
Transaxle differential brushless  
Amps:  
20 AH

Volts:  
48 Volts

Dimensions:  
56" x 27" x 52"

Seat folded down  
46"

Floor Clearance:  
6"

Batteries:  
(4) Sealed Lead Acid Maintenance Free

Forward/Reverse switch:  
Yes

High/ Low Speed switch:  
Adjustable Speed Knob / Rheo-Stat

Front tire size:  
14" tire type: 90/70-10 Tubeless

Rear tire size:  
16" tire type: 3.0 – 10 Tubeless

Charger:  
48-volt Smart Charger included

Throttle Type:  
Variable Speed Control – Twist Throttle

Key Start:  
Yes

Braking System:  
Front and Rear drum brakes with Electric Brake Assist

Drive System:  
700 Watt Brushless Motor

Headlight:  
12v LED Lights

Rear lights:  
12v Bulb

Carton Size:  
Shipped Fully assembled on Pallet – 63"x 36"x 54" (LxWxH)

Scooter Weight:  
257lbs with batteries

Battery Indicator:  
Yes

Front:  
Front Frame Suspension Bars

Rear:  
Rear Shock Absorbers

Reflectors:

N/A

Warranty:

Limited – 3 Year Warranty