


BILLING/SHIPPING INFORMATION				
Date:			Name:	
PO Number:			Address:	
Branch:				
Mark for:			City:	
Dealer Contact:			State:	
Email Address:			Zip Code:	
Phone:			Phone:	
QUOTE			<input type="checkbox"/>	ORDER



**The Dreama™ is a support system for effective 24hr postural management, suitable for children and adults**



**KEY FEATURES**

- Can be used on any hard surface-including standard bed frames, profiling beds and certain safety beds
- Allows ventilation and airflow
- Easy to clean and hygienic
- Can be used to position the user in an abducted lying position to aid in hip development
- Infinitely adjustable due to the uniquely designed glide lock system



**Information on the Jenx Dreama™ and Sleepsafe® safety bed compatibility**

The twin-sized Sleepsafe bed becomes a suitable foundation for the Dreama positioning system when certain modifications are made by Sleepsafe.

Is the Dreama intended to be used with a Sleepsafe bed?	<b>YES</b>	<b>NO</b>
	<input type="checkbox"/>	<input type="checkbox"/>

Fitting a Sleepsafe bed with a Dreama postural support system requires 2 quotes

- 1) one quote from Adaptive Imports for the Dreama
- 2) one quote from Sleepsafe which would include these modifications and a customized bed deck.

Sleepsafe quote number	
------------------------	--

Quote numbers that include ATP and patient reference must be provided to each vendor at time of quoting to ensure proper and timely collaboration.

	Twin Size
Age range (approx)	9 mo - adult
Footprint (Including cushions)	79" x 35"
Height overall	7"
Single cushion size	35" x 4" x 6"
Double cushion size	35" x 8" x 6"
Max user weight	242 lbs

**Suggested HCPCS Code - E1399\***

*\*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

## REQUIRED COMPONENTS

<b>Dreama™ Postural Support System</b>		<p style="text-align: right;"><b>Twin size includes:</b></p> <ul style="list-style-type: none"> <li>flexible aluminum base</li> <li>16 single and 2 double support cushions covered with a heat sealed, fire retardant, anti-bacterial fabric</li> <li>set of terrycloth covers for cushions- plus spares</li> <li>pair of anchor straps</li> </ul>			
		<table border="1" style="float: right; border-collapse: collapse;"> <tr> <td style="padding: 5px;"><b>DC15</b></td> <td style="padding: 5px;"><b>\$7,282</b></td> <td style="text-align: center; width: 30px;"><input type="checkbox"/></td> </tr> </table>	<b>DC15</b>	<b>\$7,282</b>	<input type="checkbox"/>
<b>DC15</b>	<b>\$7,282</b>	<input type="checkbox"/>			

## PRE-CONFIGURED KITS

Pre-configured kits allow for easy prescription of the Dreama™ postural system. Full details of support pad sizes can be found on the custom configuration page. You may wish to include additional support pads and other positioning items - these can be identified on the custom configuration page.

<p style="font-size: small; text-align: center;">*Bed frame not included or available</p>	<p><b>Side Lying Kit:</b></p> <p>Includes:</p> <ul style="list-style-type: none"> <li>1 x DM08 - "A" <b>Small positioning pad</b> with 1 x DU22 cover</li> <li>2 x DM09 "B" <b>Medium positioning pads</b> with 2 x DU25 covers</li> <li>2 x DM10 "C" <b>Large positioning pads</b> with 2 x DU28 covers</li> <li>1 x DM11 "D" <b>Extra Large positioning pad</b> with 1 x DU31 cover</li> <li>1 x DP06 "K" <b>Size 2 side lying pad</b> with 1 x DU132 cover</li> <li>1 x DU141 "M" <b>Side lying strap</b> with 1 x DU143 cover</li> </ul>				
	<table border="1" style="float: right; border-collapse: collapse;"> <tr> <td style="padding: 5px;"><b>Side Lying Kit</b></td> <td style="padding: 5px;"><b>SL</b></td> <td style="padding: 5px;"><b>\$3,915</b></td> <td style="text-align: center; width: 30px;"><input type="checkbox"/></td> </tr> </table>	<b>Side Lying Kit</b>	<b>SL</b>	<b>\$3,915</b>	<input type="checkbox"/>
<b>Side Lying Kit</b>	<b>SL</b>	<b>\$3,915</b>	<input type="checkbox"/>		
	<p><b>Supine Kit 1:</b></p> <p>Includes:</p> <ul style="list-style-type: none"> <li>2 x DM08 - "A" <b>Small positioning pad</b> with 2 x DU22 covers</li> <li>4 x DM10 "C" <b>Large positioning pad</b> with 4 x DU28 covers</li> <li>1 x DL06 "F" <b>Abduction block size 1</b> with 1 x DU52 cover</li> <li>1 x DB03 "N" <b>Large hoop</b> with 1 x DU16 cover</li> </ul>				
	<table border="1" style="float: right; border-collapse: collapse;"> <tr> <td style="padding: 5px;"><b>Supine Kit 1</b></td> <td style="padding: 5px;"><b>SP1</b></td> <td style="padding: 5px;"><b>\$3,685</b></td> <td style="text-align: center; width: 30px;"><input type="checkbox"/></td> </tr> </table>	<b>Supine Kit 1</b>	<b>SP1</b>	<b>\$3,685</b>	<input type="checkbox"/>
<b>Supine Kit 1</b>	<b>SP1</b>	<b>\$3,685</b>	<input type="checkbox"/>		
	<p><b>Supine Kit 2:</b></p> <p>Includes:</p> <ul style="list-style-type: none"> <li>2 x DM09 "B" <b>Medium positioning pad</b> with 2 x DU25 covers</li> <li>1 x DL03 "H" <b>Height adjustable T-roll</b> with 1 x DU61 cover</li> </ul>				
	<table border="1" style="float: right; border-collapse: collapse;"> <tr> <td style="padding: 5px;"><b>Supine Kit 2</b></td> <td style="padding: 5px;"><b>SP2</b></td> <td style="padding: 5px;"><b>\$2,157</b></td> <td style="text-align: center; width: 30px;"><input type="checkbox"/></td> </tr> </table>	<b>Supine Kit 2</b>	<b>SP2</b>	<b>\$2,157</b>	<input type="checkbox"/>
<b>Supine Kit 2</b>	<b>SP2</b>	<b>\$2,157</b>	<input type="checkbox"/>		

Full set of crib sized covers	<b>DU06</b>	<b>\$399</b>	<b>QTY</b>	<input type="checkbox"/>
Full set of twin sized covers	<b>DU03</b>	<b>\$950</b>	<b>QTY</b>	<input type="checkbox"/>
Additional double section for DC15	<b>DM04</b>	<b>\$324</b>	<b>QTY</b>	<input type="checkbox"/>
Individual cover for DM04	<b>DU43</b>	<b>\$121</b>	<b>QTY</b>	<input type="checkbox"/>
Additional Single Section for DC15	<b>DM03</b>	<b>\$253</b>	<b>QTY</b>	<input type="checkbox"/>
Individual Cover for DM03	<b>DU40</b>	<b>\$83</b>	<b>QTY</b>	<input type="checkbox"/>

Customized kits available upon request. Please contact Adaptive Imports.

**CUSTOM CONFIGURATION**

CUSTOM CONFIGURATION													
Slim Positioning Pad (E)		<b>Slim support pad</b> Ideal for supporting the trunk allowing for more comfortable arm positioning	Optional <b>terrycloth cover</b>	Optional terrycloth cover- padded	<b>KEY</b> + Intended use = Possible to use dependent on needs and child size x Not recommended	Supine Lying	Trunk	Pelvis	Knees-Legs	Slim Positioning Pad			
							+	=	=				
							Side Lying	Trunk	Pelvis		Knees /Legs		
							x	x	x				
DM36	DU96	DU95											
\$502	\$84	\$243											
A: 6" B: 4" C: .8"	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>							
Small Positioning Pad (A)		<b>Small positioning pad</b> Made from hygienic PU and support head rotates incrementally	Optional <b>terrycloth cover</b>	<b>Memory foam pad</b> Use in <b>addition</b> with DM08 to provide a more padded contact surface.	<b>Terrycloth cover</b> when using DM08 + DP01	Supine Lying	Trunk	Pelvis	Knees /Legs	Small Positioning Pad			
							+	+	=				
							Side Lying	Trunk	Pelvis		Knees /Legs		
							=	=	x				
DM08	DU22	DP01	DU70										
\$315	\$77	\$206	\$77										
A: 5" B: 4" C: 4" D: 4" E: 1"	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>					
Medium Positioning Pad (B)		<b>Medium positioning pad</b> Made from hygienic PU and support head rotates incrementally	Optional <b>terrycloth cover</b>	<b>Memory foam pad.</b> Use in <b>addition</b> with DM09 to provide a more padded contact surface.	<b>Terrycloth cover</b> when using DM09 + DP02	Supine Lying	Trunk	Pelvis	Knees /Legs	Medium Positioning Pad			
							=	+	+				
							Side Lying	Trunk	Pelvis		Knees /Legs		
							+	+	=				
DM09	DU25	DP02	DU73										
\$361	\$77	\$210	\$78										
A: 6" B: 4" C: 4" D: 5" E: 1"	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>					
Large Positioning Pad (C)		<b>Large positioning pad</b> Made from hygienic PU and support head rotates incrementally	Optional <b>terrycloth cover</b>	<b>Memory foam pad</b> Use in <b>addition</b> with DM10 to provide a more padded contact surface.	<b>Terrycloth cover</b> when using DM10 + DP03	Supine Lying	Trunk	Pelvis	Knees /Legs	Large Positioning Pad			
							=	+	+				
							Side Lying	Trunk	Pelvis		Knees /Legs		
							+	+	+				
DM10	DU28	DP03	DU76										
\$380	\$86	\$230	\$101										
A: 7" B: 8" C: 4" D: 6" E: 1"	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>					
X-Large Positioning Pad (D)		<b>X-Large positioning pad</b> Made from hygienic PU and support head rotates incrementally	Optional <b>terrycloth cover</b>	<b>Memory foam pad.</b> Use in <b>addition</b> with DM11 to provide a more padded contact surface.	<b>Terrycloth cover</b> when using DM11 + DP04	Supine Lying	Trunk	Pelvis	Knees /Legs	X-Large Positioning Pad			
							x	=	+				
							Side Lying	Trunk	Pelvis		Knees /Legs		
							+	+	+				
DM11	DU31	DP04	DU94										
\$441	\$101	\$237	\$113										
A: 11" B: 8" C: 4" D: 10" E: 2"	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>	QTY	<input type="checkbox"/>					

**CUSTOM CONFIGURATION- CONTINUED**

<b>Abduction Block (F)</b>		<b>Abduction Block</b> Used to ensure hips are abducted and the body is straight. Support head rotates incrementally	Optional <b>Terrycloth cover</b>	<b>Abduction Block (G)</b>		<b>Abduction Block</b> Used to ensure hips are abducted and the body is straight. Support head rotates incrementally	Optional <b>Terrycloth cover</b>	
		<b>DL06</b>	<b>DU52</b>				<b>DL07</b>	<b>DU55</b>
		<b>\$570</b>	<b>\$133</b>				<b>\$602</b>	<b>\$133</b>
	A: 3" B: 8" C: 8" D: 5"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>		A: 3" B: 9" C: 8" D: 5"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>	
<b>Knee Cups (I)</b>		<b>Knee Cups</b> Comfortably holds knees in a specified position. Support head rotates incrementally. Pair	Optional <b>Terrycloth cover. Pair</b>	<b>Knee Cups (J)</b>		<b>Knee Cups</b> Comfortably holds knees in a specified position. Support head rotates incrementally. Pair	Optional <b>Terrycloth cover. Pair</b>	
		<b>DL01</b>	<b>DU34</b>				<b>DL02</b>	<b>DU37</b>
		<b>\$578</b>	<b>\$239</b>				<b>\$610</b>	<b>\$249</b>
	A: 5" B: 3" C: 4"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>		A: 5" B: 3" C: 4"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>	
<b>Height Adjustable T-Roll (H)</b>		<b>Height Adjustable T-Roll</b> Offers knee support and prevents scissoring. Adjustable with removable pommel	Optional <b>Terrycloth cover</b>	<b>Angle Adjustable Footplate (P)</b>		<b>Angle Adjustable Footplate</b> Supports feet in supine to prevent feet dropping (each). Support head rotates incrementally	Optional <b>Terrycloth cover</b>	
		<b>DL03</b>	<b>DU61</b>				<b>DL05</b>	<b>DU67</b>
		<b>\$1014</b>	<b>\$271</b>				<b>\$432</b>	<b>\$43</b>
		QTY <input type="checkbox"/>	QTY <input type="checkbox"/>		A: 11" B: 8"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>	
<b>Side Lying Pad (K)</b>		<b>Side Lying Pad</b> encourages a neutral leg separation when side-lying. Can be used in conjunction with Side Lying Strap	Optional <b>Terrycloth cover</b>	<b>Side Lying Pad (L)</b>		<b>Side Lying Pad</b> encourages a neutral leg separation when side-lying. Can be used in conjunction with Side Lying Strap	Optional <b>Terrycloth cover</b>	
		<b>DP06</b>	<b>DU132</b>				<b>DP07</b>	<b>DU133</b>
		<b>\$590</b>	<b>\$199</b>				<b>\$593</b>	<b>\$209</b>
	A: 6" B: 18" C: 2" D: 4" E: 11" F: 3"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>		A: 6" B: 22" C: 3" D: 4" E 14" F: 4"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>	
<b>Side Lying Strap (M)</b>		<b>Side Lying Strap</b> Can be used to secure the upper leg when side lying, Requires use with 2 support pads	Optional <b>Terrycloth cover</b>	<b>Locking Feet</b>		<b>Locking Feet</b> Marks position where an accessory should be placed and keeps in place for users with high movement		
		<b>DU141</b>	<b>DU143</b>				<b>DL09</b>	
		<b>\$273</b>	<b>\$110</b>				<b>\$73</b>	
	A: 16" B: 6"	QTY <input type="checkbox"/>	QTY <input type="checkbox"/>		Set of 6	QTY <input type="checkbox"/>		
<b>Large Hoop (O)</b>		<b>Large Hoop</b> Keeps bedclothes away from the body or allow movement within a limited range	Optional <b>Terrycloth cover</b>	<b>Small Hoop (N)</b>		<b>Small Hoop</b> used to keep bedclothes away from the body or allow movement within a limited range	Optional <b>Terrycloth cover</b>	
		<b>DB03</b>	<b>DU16</b>				<b>DB01</b>	<b>DU10</b>
		<b>\$282</b>	<b>\$55</b>				<b>\$274</b>	<b>\$55</b>
		QTY <input type="checkbox"/>	QTY <input type="checkbox"/>			QTY <input type="checkbox"/>	QTY <input type="checkbox"/>	