



## **EW-M45**

The EW-M45 is a power wheelchair designed with travel in mind. Weighing just 59 pounds, this chair folds in half in only a few steps. Features include back seat pocket for storage while on-the-go, a specialty foam seat design for extreme comfort, and a large foot platform. This innovative power chair gives you the freedom to travel without having to worry about vehicle lifts or disassembling a scooter.

### EW-M45 Specifications

Maximum Speed  
3.7mph

Weight Capacity  
400 lbs.

Travel Range  
15.5 miles

Unfolding Size  
38"x23"x37"

Folding Size  
23"x13"x31"

Motor Size  
Two 180W Brushless Motors

Battery Size  
Two 12V/6AH Lithium Batteries

Brake  
Intelligent Electromagnetic Brake

Wheels:  
Front 8" solid  
Rear 12.5" solid

Seat Width  
19"

Seat Depth  
17"

Seat Height  
18"

Turning Radius  
31"

Ground Clearance  
3" at Motor

Weight w/ Batteries  
59lbs.