

Lift N' Go Electric Lift

Model 210



Standard Features

- Automatically folds up to stowaway position
- Switch box allows for easy operation of the lift
- Includes 2 retractable ratchets for securing the chair
- Includes wire harness with 30 amp resettable circuit breaker
- **SWING AWAY OPTION AVAILABLE**

Product Specifications

- All steel, e-coated, & powder finish
- 350 lb weight capacity
- 48" long X 28" wide usable platform area
- Weighs 103 lbs
- Can accommodate Class II or Class III hitch

Made in the U.S.A!

WheelChair Carrier

Toll Free: (800) 541-3213

Phone: (734) 568-6084

www.wheelchaircarrier.com



Our most popular electric lift that universally accommodates scooters and power chairs.
Perfect for effortless transportation!



*Swing Away Option
allows tailgate or
trunk access*

