

# Car Seat (MPS)

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an X the items you wish to purchase.

e-mail to: [orders@specialtomato.com](mailto:orders@specialtomato.com) / call in your order to our Pediatric Helpline at (866) 529-8407.

Account #: \_\_\_\_\_ Ship to: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_  
 Attention: \_\_\_\_\_ Phone: (     ) \_\_\_\_\_

## MPS Car Seat Features:

- FMVSS213 Compliant
- LATCH System
- FAA approved for aircraft travel for seat purchased after April 1, 2006
- Weight Capacity:
  - Small: 20 - 80 lbs.
  - Large: up to 130 lbs.
- 8" - 17" Seat Depth



## Select Your Cushions: Choose One of Each

MPS Seat Specifications	Small MPS Seat Cushions			Large MPS Seat Cushions			
	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Weight	80 lbs.			130 lbs.			
Seat Width	9"	10"	10"	11"	12"	13"	13"
Seat Depth	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"			8" - 18"			
Hip Flex Angle	-15° to +15°			-15° to +15°			
Seat Tilt Angle Forward	10°			10°			
Seat Tilt Angle Back	-30°			-30°			
Floor to Seat Height	14" - 28.5"			17" - 28.5"			

MPS Back Specifications	Small Back Cushions		Large Back Cushions		
	MS2	MS4	ML2	ML4	ML6
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult	10 yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult	48" - Adult
Maximum User Weight	80 lbs.		130 lbs.		
Back Width	8"	9"	10"	11"	13"
Back Height	8.5"	11"	11"	13"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"	26" - 30"

# Car Seat (MPS)

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an X the items you wish to purchase.

e-mail to: [orders@specialtomato.com](mailto:orders@specialtomato.com) or call in your order to our Pediatric Helpline at (866) 529-8407.

## Customize Your Size



**Small Shell (MS)**

### Shells

Choose a Shell size to determine your options




**Large Shell (ML)**

### Options

#### Headrests Cushions for both Small (MS) and Large Shells (ML)

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)



7"W x 6"H

MBH (Basic)



7.5"W x 6.5"H  
1.5"D

MBL (Basic w/ Laterals)



7"W x 6"H  
2.5"D

#### Back Cushions for Small Shell (MS)



8"W x 8.5"H



9"W x 11"H



10"W x 11"H

#### Back Cushions for Large Shell (ML)



10"W x 11"H



11"W x 13"H



13"W x 13"H

#### Seat Cushions for Small Shell (MS)



7"W x 6"D



9"W x 8"D



10"W x 10"D



10"W x 12"D


#### Seat Cushions for Large Shell (ML)



11"W x 11"D



12"W x 13"D



13"W x 15"D



13"W x 17"D

# Car Seat (MPS)

Online Dealer Resources: Get Quotes Instantly & Get Immediate Order Confirmations!

Please mark clearly with an X the items you wish to purchase.

e-mail to: [orders@specialtomato.com](mailto:orders@specialtomato.com) or call in your order to our Pediatric Helpline at (866) 529-8407.

### Car Seat Shells

Size	Part #	MSRP
Small (MS) <input type="checkbox"/>	79040109	\$ 1,371.77
Large (ML) <input type="checkbox"/>	79040309	\$ 2,507.92



Small Shell (MS)



Large Shell (ML)

## Select MPS Cushions: One of Each



### Headrest Cushions (Required) Fits both sizes of MPS Car Seat

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MHS <input type="checkbox"/>	79100051	79100052	79100054	79100059	\$ 248.10
MHL <input type="checkbox"/>	79110051	79110052	79110054	79110059	\$ 309.14
MBH (Basic Low Profile) <input type="checkbox"/>	79130051	79130052	79130054	79130059	\$ 248.10
MBL (Basic Low-Profile w/ Laterals) <input type="checkbox"/>	79140051	79140052	79140054	79140059	\$ 309.14

### Back Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS2 <input type="checkbox"/>	79120551	79120552	79120554	79120559	\$ 347.83
MS4 <input type="checkbox"/>	79120151	79120152	79120154	79120159	\$ 378.52
ML2 <input type="checkbox"/>	79120751	79120752	79120754	79120759	\$ 433.37
ML4 <input type="checkbox"/>	79120351	79120352	79120354	79120359	\$ 445.00
ML6 <input type="checkbox"/>	79120451	79120452	79120454	79120459	\$ 471.70

### Seat Cushions

Size	Cherry	Lilac	Aqua	Dark Gray	MSRP
MS3 <input type="checkbox"/>	79050151	79050152	79050154	79050159	\$ 378.52
MS5 <input type="checkbox"/>	79050251	79050252	79050254	79050259	\$ 409.22
MS7 <input type="checkbox"/>	79050651	79050652	79050654	79050659	\$ 433.37
ML1 <input type="checkbox"/>	79050751	79050752	79050754	79050759	\$ 433.37
ML3 <input type="checkbox"/>	79050351	79050352	79050354	79050359	\$ 445.00
ML5 <input type="checkbox"/>	79050451	79050452	79050454	79050459	\$ 456.63
ML7 <input type="checkbox"/>	79050851	79050852	79050854	79050859	\$ 483.57

Pricing valid 1/1/2024 - 12/31/2024

\$

Your MSRP Total

Toll-Free Helpline: (866) 529-8407 e-mail: [orders@specialtomato.com](mailto:orders@specialtomato.com)

Order Forms Available Online at: [www.specialtomato.com/orderforms.html](http://www.specialtomato.com/orderforms.html)

OFCAR01-24LIT