

# 2025 order form for medium

Use dimension chart to select appropriate size.

Scan for sizing and pricing:  
<https://www.rifton.com/activitychair>



		Standard R840 medium		Hi/lo R850 medium		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$695 <input type="checkbox"/>	R851	\$1,440 <input type="checkbox"/>
		Adjustable w/ back spring	R842	\$790 <input type="checkbox"/>	R852	\$1,565 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R843	\$400 <input type="checkbox"/>	R853	\$2,950 <input type="checkbox"/>
		With spring	R844	\$480 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$280 <input type="checkbox"/>	R868	\$280 <input type="checkbox"/>
		Forearm prompts	R869	\$700 <input type="checkbox"/>	R869	\$700 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R893	\$70 <input type="checkbox"/>	N/A	
		Long	R894	\$90 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$340 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$375 <input type="checkbox"/>	N/A	
* Pads	CHOOSE	Tan	R835	\$285 <input type="checkbox"/>	R835	\$285 <input type="checkbox"/>
		Blue	R836	\$285 <input type="checkbox"/>	R836	\$285 <input type="checkbox"/>
		Red	R837	\$285 <input type="checkbox"/>	R837	\$285 <input type="checkbox"/>
		Pink	R974	\$285 <input type="checkbox"/>	R974	\$285 <input type="checkbox"/>
		Purple	R975	\$285 <input type="checkbox"/>	R975	\$285 <input type="checkbox"/>
		Green	R976	\$285 <input type="checkbox"/>	R976	\$285 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897	\$270 <input type="checkbox"/>	R897	\$270 <input type="checkbox"/>
		Flat	R898	\$270 <input type="checkbox"/>	R898	\$270 <input type="checkbox"/>
		Adjustable winged	R881	\$300 <input type="checkbox"/>	R881	\$300 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828	\$230 <input type="checkbox"/>	R828	\$230 <input type="checkbox"/>
		Small with chest strap	R829	\$285 <input type="checkbox"/>	R829	\$285 <input type="checkbox"/>
		Large	R838	\$244 <input type="checkbox"/>	R838	\$244 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)			R878	\$170 <input type="checkbox"/>	Included	
Tray (requires armrests)			R668	\$310 <input type="checkbox"/>	R668	\$310 <input type="checkbox"/>
Handhold			R649	\$32 <input type="checkbox"/>	R649	\$32 <input type="checkbox"/>
Chest strap (wide)			R863	\$105 <input type="checkbox"/>	R863	\$105 <input type="checkbox"/>
Butterfly harness	CHOOSE	Regular	R864	\$140 <input type="checkbox"/>	R864	\$140 <input type="checkbox"/>
		Slim cut	R882	\$140 <input type="checkbox"/>	R882	\$140 <input type="checkbox"/>
Thigh belt			R865	\$110 <input type="checkbox"/>	R865	\$110 <input type="checkbox"/>
Pelvic harness			R866	\$200 <input type="checkbox"/>	R866	\$200 <input type="checkbox"/>
Hip guides (pair)			R867	\$240 <input type="checkbox"/>	R867	\$240 <input type="checkbox"/>
Abductor	CHOOSE	Small	R884	\$145 <input type="checkbox"/>	R884	\$145 <input type="checkbox"/>
		Large	R891	\$165 <input type="checkbox"/>	R891	\$165 <input type="checkbox"/>
Adductors (pair)			R892	\$230 <input type="checkbox"/>	R892	\$230 <input type="checkbox"/>
Leg prompt			R858	\$255 <input type="checkbox"/>	R858	\$255 <input type="checkbox"/>
Footboard (recommended for chairs with long legs)			R848	\$545 <input type="checkbox"/>	Included	
Ankle straps (pair)			R899	\$105 <input type="checkbox"/>	R899	\$105 <input type="checkbox"/>
Sandals (pair)			R954	\$110 <input type="checkbox"/>	R954	\$110 <input type="checkbox"/>
Wedges (pair, require sandals)			R698	\$85 <input type="checkbox"/>	R698	\$85 <input type="checkbox"/>
Lumbar and seat support kit			R887	\$90 <input type="checkbox"/>	R887	\$90 <input type="checkbox"/>
Backrest filler pad			R849	\$140 <input type="checkbox"/>	R849	\$140 <input type="checkbox"/>
Whitmyer headrest adapter			R804	\$65 <input type="checkbox"/>	R804	\$65 <input type="checkbox"/>
Directional Assist Wheel retrofit kit			N/A		R965	\$290 <input type="checkbox"/>
					R966 (5 kits)	\$1,450 <input type="checkbox"/>

\* Required



# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



## Medium chair at its largest dimensions


- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




**Medium standard base**  
R840 Rifton Activity Chair



**Medium hi/lo base**  
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R840 medium Standard base
Frame width	short legs: 23¾
	long legs: 24½
	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above floor	short legs: 12½–16½
	long legs: 19–23
	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½ 14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lb)	150

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½ 14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lb)	150